


# msdi Drilling Kit

Ref: MDRSK

**DRSP15** – Spade Drill 1.5  
**DRLMK** – Marking Drill Ø1.9

**DRL40** – Drill Ø4.0 L16  
**DRL45** – Drill Ø4.5 L16  
**DRL50** – Drill Ø5.0 L16  
**DRL55** – Drill Ø5.5 L16



**DRL20** – Drill Ø2.0 L16  
**DRL25** – Drill Ø2.5 L16  
**DRL28** – Drill Ø2.8 L16  
**DRL32** – Drill Ø3.2 L16  
**DRL365** – Drill Ø3.65 L16



Please scan for more  
information for use:



# Drilling Protocol


○ Throughout entire implants length

● In cases of a hard bone drill through the cortical layer


\* All measurements are in mm

\*\* The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment.

## Alef



## Til



## Conix

