


# msdi Prosthetic kit

Ref: MSKP



**RTAQ** - Torque Wrench

**DRPL** - ratchet prosthetics long 15mm

**DMPL** - Motor Prosthetics long

**DML** - manual long 15mm

**DRPS** - Ratchet prosthetics short 10mm

**DMPS** - Motor Prosthetics short

**DMS** - Manual Short 10mm



Til System Recommended Torque	
Til System	Max. Recommended Torque – Ncm
Til Implant	45
Cover Screw	25
Healing Abutment	15
Straight Abutment	25
Angled Abutment	25
Castable Abutment	15
Castable Abutment with Titanium Base	25
Ball Attachment	25
Lock Attachement	25
Multi Unit Straight Abutment	15-20
Multi Unit Angled Abutment	15-20
Multi Unit Healing Abutment	15-20
Ti-Base	25
Multi Unit Ti-Base	15-20
Multi Unit Ti-Base	15-20

Alef System Recommended Torque	
Alef System	Max. Recommended Torque – Ncm
Alef Implant	60
Cover Screw	35
Healing Abutment	20
Straight Abutment	35
Angled Abutment	35
Castable Abutment	20
Castable Abutment with Titanium Base	35
Castable Abutment with Chrome Cobalt Base	35
Ball Attachment	35
Lock Attachement	35
Multi Unit Straight Abutment	15-20
Multi Unit Angled Abutment	15-20
Multi Unit Healing Abutment	15-20
Ti-Base	35
Multi Unit Ti-Base	15-20

Conix System Recommended Torque	
Conix System	Max. Recommended Torque – Ncm
Conix Implant	60
Cover Screw	35
Healing Abutment	20
Straight Abutment	35
Angled Abutment	35
Castable Abutment	20
Castable Abutment with Titanium Base	35
Castable Abutment with Chrome Cobalt Base	35
Ball Attachment	35
Lock Attachement	35
Multi Unit Straight Abutment	15-20
Multi Unit Angled Abutment	15-20
Multi Unit Healing Abutment	15-20
Ti-Base	35
Multi Unit Ti-Base	15-20