

msdi surgical kit

Ref: MSK

DRLMK - Marking drill Ø1.9

DRPL - Driver ratchet
prosthetics long 15mm


DMPL - Driver Motor
Prosthetics long

DML - Driver manual long
15mm

DRIML - Driver ratchet
implants long 15mm

DMIL - Drivers Motor
implants long

DRLEX - Drill Extender



RTAQ - Torque Wrench

DRL20 - Drill Ø2.0 L16

DRL28 - Drill Ø2.8 L16

DRL32 - Drill Ø3.2 L16

DRL365 - Drill Ø3.65 L16

DRL40 - Drill Ø4.0L16

DRL45 - Drill Ø4.5 L16



Please scan for more
information for use



msdi Drilling Protocol

○ Throughout entire implants length

● In cases of a hard bone drill through the cortical layer

* All measurements are in mm

** The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment

